Features	Infinity Ice Effect I	Infinity Ice Effect II	Infinity Sapphire Ice Effect
	# 42 # D D D D D D D D D D D D D D D D D D		
Pulse Light Repetition Rate	Up to 1.1 Hz (1 Pulse Every 0.9 Seconds)	Up to 1.1 Hz (1 Pulse Every 0.9 Seconds)	Up to 1.1 Hz (1 Pulse Every 0.7 Seconds)
Features	Two functional modes: Hair Removal and Skin Rejuvenation	Hair removal with integrated skin rejuvenation functionality — no need to switch modes	Hair removal with three independent sapphire filters and integrated skin rejuvenation functionality — no need to switch modes
Skin Rejuvenation Wavelength	530 – 750 nm	530 – 750 nm	470 – 1200 nm
Hair Removal Wavelength	550 – 1200 nm	550 – 1200 nm	Face: 550nm – 1200nm. Body: 470nm – 1200nm. Bikini Area (VIO): 600nm – 1200nm.
Lamp Window Size	1.6 inches – 4 cm ²	1.6 inches – 4 cm ²	1.259 inches – 3.2 cm ²

Maximum Power	19 J/cm ²	23 J/cm ²	24 J/cm ²
Adjustable Energy Levels	5 Levels of Adjustable Energy	10 Levels of Adjustable Energy	3 Levels of Adjustable Energy
Precision Unit Light Window Size	X	0.79 inches – 2 cm ²	X
Cooling Effect	Ice Effect	Ice Effect	Ice Effect
Universal Voltage	110-220V AC~50/60Hz	110-220V AC~50/60Hz	110-220V AC~50/60Hz
Power Consumption	<43W	<43W	<43W

Power Cord	✓	✓	√
Suitable For	Unisex and suitable for all types of hair colors and skin types	Unisex and suitable for all types of hair colors and skin types	Unisex and suitable for all types of hair colors and skin types
Dual Treatment	Yes, for permanent hair removal and skin rejuvenation	Yes, for permanent hair removal and skin rejuvenation	Yes, for permanent hair removal and skin rejuvenation

Please note: The final effectiveness, regardless of the treatment duration with one model or another, is the same across all three products.

Key features and specifications shared by the Infinity Ice Effect I, II and Sapphire

- 1. Universal Voltage: 110-220V AC~50/60Hz.
- 2. Power consumption: <43W.
- 3. It works with a power cord.
- 4. Cooling Effect: Ice Effect.
- 5. Unisex and suitable for all types of hair colors and skin types.
- 6. Effectiveness: The final effectiveness, regardless of the treatment duration with one model or another, is the same across all three products.
- 7. Dual treatment, for permanent hair removal and skin rejuvenation.