

LESCOLTON® INFINITY ICE EFFECT II – POWERED BY ELOS ME®
USER GUIDE

Elos Me®

LESCOLTON® INFINITY ICE EFFECT II T012C – POWERED BY ELOS ME®
USER GUIDE



1800 North Shelby Oaks Drive
Memphis, TN 38134

Questions? Get in touch with us via live chat or email us:
hello@elosme.com

CONTENTS

Product Introduction	03
Contraindications	03
Considerations	03
Safety warnings	03
Safety operating regulations	03
How to use	03
Treatment intervals	04
Specifications	05
Table of simple troubleshooting	05
Contact Information	05

INTRODUCTION

- The Lescolton Infinity Icecool II powered by Elos Me is a dual system with the functions of hair removal and skin rejuvenation. It integrates the multi-term application of optical e-light technology in the field of Beauty & Health.
 - Before using this device, please carefully read the instructions.
 - Suitable for both men and women.
 - Individuals with dark skin should begin with the lowest level and then adjust to the appropriate level based on the adaptability of their skin.
-

CONTRAINDICATIONS

- Not suitable for individuals under 16 years old.
 - Avoid use in cases of skin inflammation.
 - Avoid use in cases of photosensitivity.
 - Avoid use if you have an implanted medical electronic device that may lead to operational errors.
-

CONSIDERATIONS

- This product has not been tested on pregnant women.
 - Do not use the device directly on tattoos, as they may absorb excessive energy.
 - Avoid use around the eyes, neck, or nipples, as the skin in these areas is particularly sensitive.
 - If the skin remains inflamed or swollen after 24 hours, discontinue use until symptoms subside.
 - Do not use the device directly on moles or skin tags. It is recommended to cover them with a white sticker or cosmetic pencil before treatment to prevent irritation or potential skin reactions.
-

SAFETY WARNINGS

- Do not obstruct the vents during use.
 - Ensure that the device does not come into contact with water or moisture.
 - It is prohibited to use the device in flammable or explosive conditions.
-

SAFETY OPERATING REGULATIONS

- Avoid aiming the light at your eyes.
 - Do not look directly at the light.
 - Before using the device, ensure that you have removed all metallic objects such as rings and bracelets to prevent device malfunction.
 - Use the machine in a room with a temperature below 25 degrees Celsius.
-

HOW TO USE

- Shave or wax the hair around the treatment area before using the device.
- Insert the power cable into the machine and plug it in.
- Turn on the machine using the touch power button on the left side of the LCD screen.

Press it for 2 or 3 seconds until the device turns on.

- Use the same power button to select the desired power level. Each time you press the button, the power level will increase. If you do not feel any tingling or discomfort at the lowest power levels, it's best to select the highest possible power level.
- Place the lamp on your skin, and when the green light on the left side of the device flashes, press the fire button on the handle of the device. You can repeat shots in the same area 2 or 3 times for better results. To issue the next shot, wait for the green light on the left side of the machine to flash again.
- To enable the machine to shoot automatically without pressing the shutter button, you can enter the automatic use mode. To enter the automatic mode, press the shutter button for two or three seconds until you hear two beeps. This way, the machine will emit shots automatically when the lamp is resting on the skin. Keep the lamp supported until 2 or 3 shots are emitted for the best results.
- To return to manual shooting mode, press the shutter button again for 2 to 3 seconds until you hear two beeps.
- To use the facial adapter, simply attach the adapter to the machine's lamp.



TREATMENT INTERVALS

- During the first month, treat every 3 to 5 days.
- In the second month, treat weekly.
- In the third month, treat every 15 days.
- Starting from the fourth month, treat every 30 days until permanent hair removal is achieved.

Permanent hair removal is a process that takes several months. Immediate results are not expected for two reasons:

1. All body hair must be in the active growth stage to be affected by the treatment, which may take several months.
2. Multiple treatments are necessary to weaken the hair root in the active growth stage until it is permanently affected, preventing hair regrowth. Visible changes will occur over time.

Individuals with blond, grey, or red hair typically require twice as many sessions as those with black or brown hair. Therefore, the treatment duration may extend to 6 to 12 months.

SPECIFICATIONS

- Main Unit Size: 182x78x151mm
- Package (Box) Size: 231x191x93mm
- Power Supply: Input: 100-240V AC-50/60Hz
- Weight: 1.1 kg

- Irradiation Area: 4 cm²
- External Environment: Temperature Range: -20°C to 50°C
- Warranty: 6-year warranty

TABLE OF SIMPLE TROUBLESHOOTING

- No Power: Ensure that the power cord is securely connected.
- Can't Turn On the Machine: Press and hold the "on" button for 2 seconds.
- Level Indicator Is Flashing After Boot: The cartridge may not be properly installed.
- Automatic Shutdown: The machine automatically shuts down after 5 minutes of inactivity.

Please follow all safety instructions in this manual to ensure safe and effective use of your device. For any questions or concerns, contact our customer support.

CONTACT INFORMATION

Elos Me & Lescolton
1800 North Shelby Oaks Drive
Memphis, TN 38134
Email: hello@elosme.com
Website: www.elosme.com

This product complies with international safety standards.

© Elos Me. All rights reserved. No part of this user guide may be reproduced or transmitted in any form without the prior written permission of Elos Me.