

**LESCOLTON® INFINITY ICE EFFECT II – POWERED BY ELOS ME®**  
**USER GUIDE**

**Elos Me®**

**LESCOLTON® INFINITY ICE EFFECT II T012C – POWERED BY ELOS ME®**  
**USER GUIDE**



1800 North Shelby Oaks Drive  
Memphis, TN 38134

Questions? Get in touch with us via live chat or email us:  
[hello@elosme.com](mailto:hello@elosme.com)

**CONTENTS**

Product Introduction .....03  
Contraindications .....03  
Considerations .....03  
Safety warnings .....03  
Safety operating regulations .....03  
How to use .....03  
Treatment intervals .....04  
Specifications .....05  
Table of simple troubleshooting .....05  
Contact Information .....05

## **INTRODUCTION**

- The Lescolton Infinity Icecool II powered by Elos Me is a dual system with the functions of hair removal and skin rejuvenation. It integrates the multi-term application of optical e-light technology in the field of Beauty & Health.
  - Before using this device, please carefully read the instructions.
  - Suitable for both men and women.
  - Individuals with dark skin should begin with the lowest level and then adjust to the appropriate level based on the adaptability of their skin.
- 

## **CONTRAINDICATIONS**

- Not suitable for individuals under 16 years old.
  - Avoid use on inflamed skin.
  - Avoid use in cases of photosensitivity.
  - Avoid use if you have an implanted medical electronic device that may lead to operational errors.
- 

## **CONSIDERATIONS**

- This product has not been tested on pregnant women.
  - Do not use the product directly on tattoos, as they may absorb a significant amount of energy.
  - Avoid use around the eyes, neck, or nipples, as the skin in these areas is relatively sensitive.
  - If the skin continues to be inflamed or swollen after 24 hours, discontinue use until symptoms disappear.
- 

## **SAFETY WARNINGS**

- Do not obstruct the vents during use.
  - Ensure that the device does not come into contact with water or moisture.
  - It is prohibited to use the device in flammable or explosive conditions.
- 

## **SAFETY OPERATING REGULATIONS**

- Avoid aiming the light at your eyes.
  - Do not look directly at the light.
  - Before using the device, ensure that you have removed all metallic objects such as rings and bracelets to prevent device malfunction.
  - Use the machine in a room with a temperature below 25 degrees Celsius.
- 

## **HOW TO USE**

- Shave or wax the hair around the treatment area before using the device.
- Insert the power cable into the machine and plug it in.

- Turn on the machine using the touch power button on the left side of the LCD screen. Press it for 2 or 3 seconds until the device turns on.
- Use the same power button to select the desired power level. Each time you press the button, the power level will increase. If you do not feel any tingling or discomfort at the lowest power levels, it's best to select the highest possible power level.
- You can choose from four treatment modes using the button on the right side of the LCD screen (work mode switch button): hair removal, hair removal with cooling, skin rejuvenation, and skin rejuvenation with cooling.
- To switch between hair removal and skin rejuvenation modes, press the "work mode switch button." HR indicates the hair removal mode, and SR indicates the skin rejuvenation mode.
- To turn on/off the Ice-cool care mode, press the "work mode switch button" and hold it down for 2/3 seconds. The snowflake icon indicates that the Ice-cool care mode is on; otherwise, it's off.
- Place the lamp on your skin, and when the green light on the left side of the device flashes, press the fire button on the handle of the device. You can repeat shots in the same area 2 or 3 times for better results. To issue the next shot, wait for the green light on the left side of the machine to flash again.
- To enable the machine to shoot automatically without pressing the shutter button, you can enter the automatic use mode. To enter the automatic mode, press the shutter button for two or three seconds until you hear two beeps. This way, the machine will emit shots automatically when the lamp is resting on the skin. Keep the lamp supported until 2 or 3 shots are emitted for the best results.
- To return to manual shooting mode, press the shutter button again for 2 to 3 seconds until you hear two beeps.
- To use the facial adapter, simply attach the adapter to the machine's lamp.




---

### **TREATMENT INTERVALS**

- During the first month, treat every 3 to 5 days.
- In the second month, treat weekly.
- In the third month, treat every 15 days.
- Starting from the fourth month, treat every 30 days until permanent hair removal is achieved.

Permanent hair removal is a process that takes several months. Immediate results are not expected for two reasons:

1. All body hair must be in the active growth stage to be affected by the treatment, which may take several months.
2. Multiple treatments are necessary to weaken the hair root in the active growth stage until it is permanently affected, preventing hair regrowth. Visible changes will occur over time.

Individuals with blond, grey, or red hair typically require twice as many sessions as those with black or brown hair. Therefore, the treatment duration may extend to 6 to 12 months.

---

### **SPECIFICATIONS**

- Main Unit Size: 182x78x151mm
  - Package (Box) Size: 231x191x93mm
  - Power Supply: Input: 100-240V AC-50/60Hz
  - Weight: 1.1 kg
  - Irradiation Area: 4 cm<sup>2</sup>
  - External Environment: Temperature Range: -20°C to 50°C
  - Warranty: 6-year warranty
- 

### **TABLE OF SIMPLE TROUBLESHOOTING**

- No Power: Ensure that the power cord is securely connected.
  - Can't Turn On the Machine: Press and hold the "on" button for 2 seconds.
  - Level Indicator Is Flashing After Boot: The cartridge may not be properly installed.
  - Automatic Shutdown: The machine automatically shuts down after 5 minutes of inactivity.
- 

Please follow all safety instructions in this manual to ensure safe and effective use of your device. For any questions or concerns, contact our customer support.

---

### **CONTACT INFORMATION**

Elos Me & Lescolton  
1800 North Shelby Oaks Drive  
Memphis, TN 38134  
Email: [hello@elosme.com](mailto:hello@elosme.com)  
Website: [www.elosme.com](http://www.elosme.com)

---

This product complies with international safety standards.

© Elos Me. All rights reserved. No part of this user guide may be reproduced or transmitted in any form without the prior written permission of Elos Me.