

**LESCOLTON® SAPPHIRE ICE EFFECT - POWERED BY ELOS ME®**  
**USER GUIDE**

**Elos Me®**

**LESCOLTON® SAPPHIRE ICE EFFECT - POWERED BY ELOS ME®  
USER GUIDE**



1800 North Shelby Oaks Drive  
Memphis, TN 38134

Questions? Get in touch with us via live chat or email at: [hello@eloame.com](mailto:hello@eloame.com)

## CONTENTS

- Product Introduction.....03
- Contraindications.....03
- Considerations.....03
- Safety Warnings.....03
- Safety Operating Regulations.....03
- How to Use the Device.....04
- Treatment Intervals.....04
- Specifications.....05
- Table of Simple Troubleshooting.....05

## **INTRODUCTION**

- The Professional Sapphire technology for home use revolutionizes the process by introducing an innovative linear ice-sensing system with three distinct sapphire filters tailored for body, bikini area, and facial hair removal.
  - Before use, please carefully read the instructions.
  - Suitable for both men and women.
  - People with dark skin should start with the lowest level and then adjust it based on their skin's adaptability.
- 

## **CONTRAINDICATIONS**

- Not suitable for individuals under 16 years old.
  - Avoid use in cases of skin inflammation.
  - Avoid use in cases of photosensitivity.
  - Do not use if you have implanted medical electronic devices that may lead to operational errors.
- 

## **CONSIDERATIONS**

- This product has not been tested on pregnant women.
  - Avoid using the product directly on tattoos, as they may absorb excessive energy.
  - Do not use around the eyes, neck, or nipples, as the skin in these areas is relatively sensitive.
  - If the skin remains inflamed or swollen after 24 hours, discontinue use until symptoms disappear.
- 

## **SAFETY WARNINGS**

- Do not obstruct the vents during use.
  - Keep the device away from moisture and water.
  - Use of the device in flammable or explosive conditions is prohibited.
- 

## **SAFETY OPERATING REGULATIONS**

- Avoid directing the light at your eyes.
- Do not look directly at the light.
- Before use, ensure that you have removed all metallic objects, such as rings or bracelets, to prevent device malfunctions.

- Use the machine in a room with a temperature below 25 degrees Celsius.
- 

## HOW TO USE

- Shave or wax the hair around the treatment area before using the device.
  - Insert the power cable into the machine and plug it in.
  - Turn on the machine using the power button on the right side of the unit by pressing it for 2 or 3 seconds until the device turns on.
  - Using the same power button, select your desired power level. Each press increases the power level. If you don't feel any tingling or discomfort at the lowest power levels, it's best to select the highest level.
  - Three treatment modes correspond to the three distinct sapphire filters for body, bikini area, and facial hair removal. Place the filter for your desired treatment (Body, Face, or VIO) on the top left-hand side of the machine.
  - Position the unit vertically with the lamp resting on the skin to prepare for light pulse emission.
  - Press the light pulse emission button called 'Shot' located on the front of the unit. For best results, repeat shots on the same area 2 or 3 times.
  - To enable automatic shooting without pressing the shutter button, you can enter automatic use mode. To do so, press the shutter button for two or three seconds until you hear two beeps. In this mode, the machine will emit shots automatically when the lamp rests on the skin. Keep the lamp supported until 2 or 3 shots are emitted for the best results.
  - To return to manual shooting mode, press the shutter button again for 2 to 3 seconds until you hear two beeps.
- 

## TREATMENT INTERVALS

Recommended treatment schedule for optimal results:

- **First month:** Every 3 to 5 days
- **Second month:** Weekly
- **Third month:** Every 15 days
- **From the fourth month onward:** Every 30 days, until you achieve permanent hair removal

Achieving permanent hair removal is a process that takes several months. Results are not immediate due to two reasons:

1. Waiting for all body hair to be in the active growth stage, allowing it to be affected by the treatment, which may take several months.

2. Several treatments are required for hair in the active growth stage to weaken the hair root permanently. Changes will become noticeable over time. People with blond, grey, or red hair typically require twice as many sessions as people with black or brown hair, which can extend the treatment duration to 6 to 12 months.

---

## **SPECIFICATIONS**

- Power Supply: Input: 100-240V AC-50/60Hz.
- Weight: 1.76 lb.
- Irradiation Area: 3.2 cm<sup>2</sup>.
- Pulse Light Repetition Rate: up to 1.43 Hz (1 pulse every 0.7 seconds).
- Professional Technology with Three Independent Sapphire Filters for Treatment Optimization, Based on the Targeted Body Area: Body, Face, and Bikini. Wavelength Light Spectrum of Each Filter: Face: 550nm – 1200nm, Body: 470nm – 1200nm, Bikini Area (VIO): 600nm – 1200nm.
- Lamp Window Size: 1.259 inches – 3.2 cm<sup>2</sup>.
- Maximum Power: 24 J/cm<sup>2</sup>.
- Energy Levels: 3 adjustable energy levels.
- External Environment: Temperature: -20 to 50 degrees Celsius.
- Warranty: 6-year warranty.
- Model: LS-T108.

---

## **TABLE OF SIMPLE TROUBLESHOOTING**

- No power: Ensure the power cord is securely connected.
- Unable to turn on the machine: Press the “on” button for 2 seconds.
- Level indicator flashing after boot: The filter may not be installed correctly.
- Automatic shutdown: The machine automatically shuts down after 5 minutes of inactivity.

---

Please follow all safety instructions in this manual to ensure safe and effective use of your device. For any questions or concerns, contact our customer support.

---

## **CONTACT INFORMATION**

Elos Me & Lescolton  
1800 North Shelby Oaks Drive

Memphis, TN 38134

Email: [hello@elosme.com](mailto:hello@elosme.com)

Website: [www.elosme.com](http://www.elosme.com)

---

This product complies with international safety standards.

---

© Elos Me. All rights reserved. No part of this user guide may be reproduced or transmitted in any form without the prior written permission of Elos Me.